

Medical, Laboratory, Pharmacy, & Life Science Cold Storage

accucold 



Purpose-Built Equipment For Institutional Applications

ACCUCOLD DIVISION, FELIX STORCH, INC. • ISO 9001:2015 CERTIFIED • 888.4.MEDLAB • INFO@ACCUCOLD.COM • ACCUCOLD.COM

accucold

The Accucold Difference

Accucold, a division of Felix Storch, Inc., specializes in manufacturing purpose-built refrigeration products for the medical, pharmaceutical, laboratory, and institutional industries. Our extensive product selection offers a range of cooling solutions, from basic locking compact refrigerators to built-in vaccine storage units specially designed to comply with **CDC and VFC guidelines**. We are **ISO 9001:2015 certified** and maintain an ISO/IEC 17025:2015 certified calibration laboratory for our monitoring equipment.

MED/LAB/PHARMACY REFRIGERATION	3
Performance Series Refrigeration Highlights	4
New Options for Optimized Storage & Security	5
Pharma-Vac Performance Series Refrigerators	6-7
Pharma-Lab Performance Series Refrigerators & Freezers	8-9
Med-Lab Performance Series Refrigerators	10-11
MED2 Series Medical All-Freezers	12-13
ACR Series Pharmaceutical Refrigerators	14
Stacked Refrigerator-Freezer Combinations.....	15
PATIENT CARE/HOSPITAL/BREAK ROOM	16
PURETHERM Blanket Warmers	17
JCAHO Compliant Nutrition Centers.....	18
Break Room Refrigerator-Freezers.....	19
MOMCUBE Series Breast Milk Refrigeration.....	20-21
Microwaves.....	22
Icemakers & Dispensers.....	23
Portable Refrigerator/Freezers	24
GENERAL PURPOSE/NSF-7 Refrigeration	25
Healthcare Series Refrigerators & Freezers	26
PLUS2 Series Refrigerators & Freezers	27
Low Temperature Freezers	28
Commercial Chest Refrigerators & Freezers	29
Options & Accessories	30-31



MED/LAB/PHARMACY REFRIGERATION

Accucold's Performance Series features purpose-built equipment specially designed and constructed to meet the demanding storage needs of medical, laboratory, vaccine, and life science applications.

Performance Series Refrigeration



Solid & glass door choices from 1 to 49 cu.ft.

Purpose-built with advanced features to help support CDC and VFC guidelines for proper vaccination storage:



Audible alarms for hi/lo temp., power out, & open door

- Keyed lock for secure interior
- Intelligent microprocessor control panel with password protection
- Optimized forced air cooling for fast temp. recovery
- Digital display of the current & hi/lo temp.
- Silver ion-coated handle provides protection from germs & microbes
- LED lighting for an efficient display
- Self-closing user-reversible door
- Adjustable plastic-coated wire shelves
- Durable interior construction from powder-coated steel



EXTENDED WARRANTY
 • 2-year parts & labor
 • 5-year compressor



ECO-FRIENDLY
 Environmentally friendly hydrocarbon refrigerants



COMPATIBLE
 with Pyxis®, Omnicell®, & AcuDose RX®



ADA OPTION
 32" high choices for use under lower ADA compliant counters

Adjustable operating control range: 2 to 8°C • Temperature variation: +/- 1°C*

*The published temperature variation is derived from the maximum deviation of an NTC sensor in a 1 oz. vial located nearest the chamber geometric center during a 24-hour test period

New Options for Optimized Storage & Security

Ventilated Drawers

Full-extension removable drawers made from powder-coated aluminum, factory-installed in any 24" wide Pharma-Vac or Med-Lab refrigerator



Suffix "DR"

Locker Compartments

Interior lockers secured by combination locks for separated storage. Staff can set their own security code and override with provided keys



Suffix "Locker"

Keyless Access Locks

No keys required with a keypad security system that allows you to program up to 19 codes to accommodate your staff. Available as a factory-installed option on all Performance Series units



Add Option "KAL"

Digital Data Loggers

NIST calibrated loggers with internal storage and removable USB drives for virtually unlimited data storage. Loggers are calibrated in our ISO/IEC 17025:2015 certified laboratory and display the current and min/max temperatures in Celsius and Fahrenheit



Add Option "DL2B"

All units can be combined with Accucold freezers using our stacking racks

*Specifications are subject to change without notice. For the most current information, visit accucold.com

Pharma-Vac Performance Series

Purpose-built for pharmacy, medication, and vaccination applications to support meeting CDC/VFC vaccine storage guidelines



Countertop



ADA



Undercounter



Full Sized



Advanced Temperature Control

- Intelligent microprocessor** digital temperature controller
- Improved temperature accuracy** with a buffered probe in glycol-filled bottle to simulate temperature of stored product
- Additional sensor** reads interior air temperature for superior temperature control
- Adjustable operating control range** from +2°C to +8°C (+35.6°F to +46.4°F)
- Digital display** of the min/max temperature (updated in 15 minute increments)
- Temperature viewable in Celsius and Fahrenheit** to the nearest tenth of a degree
- Password protection** of control parameters beyond setpoint
- Optimized forced air cooling** for excellent temperature stability and uniformity with rapid recovery after door openings and product loading



Multiple Alarms For Safer Monitoring

- Audible/visual alarms** for high/low temperature excursions
- Open door alarm** sounds if the door is left ajar (default 4 minutes with selectable 1-99 minute range)
- Audible power failure** and sensor failure alarms
- Mute function** for alarms
- Remote alarm** contacts in rear of the unit




Fully Featured Interior For Convenient Storage

- Energy efficient LED lighting** with internal on/off switch
- Adjustable plastic-coated wire shelves** for improved circulation
- Multiple shelf positions** spaced at 1/2" to allow flexible adjusting to support smaller items
- Adjustable cycle defrost** operation
- 1/2" **probe port** included to accept user-provided monitoring equipment



Quality Design For Reliable Performance

- Hospital grade cord** with a green dot plug
- Antimicrobial sliver ion-coated handle** for added protection from germs
- Self-closing door** offers added storage protection
- User-reversible door swing** (self-closing function compatible with RHD position only)
- Environmentally-friendly** natural hydrocarbon refrigerant for reduced global warming potential
- Meets EPA and SNAP** requirements
- Convenient security** with a top-mounted door lock (2 keys included)
- Magnetic gaskets** to ensure a positive seal
- Compatible with** Pyxis®, Omnicell® and AcuDose RX®



SKU: ARG12PV



Remote contacts



Temperature sensor in glycol bottle



Microprocessor digital temperature control



New Options Available

- Keyless Locks
- Interior Lockers
- Ventilated Drawers

See Page 5



ARS1PV

ARG3PV

ARG8PVDR

ARS12PV

ARG15PVLOCKER

Technical Specifications (All Units):

Defrost type: Cycle
Thermostat: Intelligent microprocessor digital
Temperature sensors: 2 (air probe & buffered probe)
Adjustable operating control range: 2 to 8°C (35.6 to 46.4°F)
Temperature variation: +/- 1°C*
Shelving: Adjustable white plastic coated wire
Voltage: 115V/60 Hz
Cord: NEMA 5-15P
Lock: Factory-installed key lock
Refrigerant: Eco-friendly R600a
Light: LED
Access port: ½"

Door: Reversible
Fan: Internal
Certifications: ETL, ETL-C
Parts & Labor Warranty: 2 years
Compressor warranty: 5 years

Included with units:
 2 keys
 Do Not Unplug removable sticker
 Do Not Adjust Temperature removable sticker
 Do Not Store Personal Items removable sticker
 User manual

DESIGNED TO PROVIDE +/-1°C VARIATION
 The published temperature variation is derived from the maximum deviation of an NTC sensor in a 1 oz. vial located nearest the chamber geometric center during a 24 hour test period.

100% GUARANTEE
Extended Warranty
 2-year parts & labor
 5-year compressor

Eco-Friendly
 Uses environmentally friendly hydrocarbons for a greener footprint

Compatible
 with Pyxis®, Omnicell®, & AcuDose RX®

Logging Option
 Add our D2LB data logger for additional assurance

	MODEL	CAPACITY	HEIGHT x WIDTH x DEPTH	INSTALLATION	DOOR	CABINET	SHELVES	AMPERAGE	SHIPPING WT.
🔒	ARS1PV	1 cu.ft.	21.5" x 17.5" x 19.75"	Freestanding	White	White	3	0.6	70 lbs
🔒	ARG1PV	1 cu.ft.	21.5" x 17.5" x 19.75"	Freestanding	Glass	White	3	0.6	73 lbs
🔒	ARS32PVBIADA	2.83 cu.ft.	32" x 18.5" x 22.63"	Recessed	White	White	4	0.6	113 lbs
🔒	ARG31PVBIADA	2.83 cu.ft.	32" x 18.5" x 22.63"	Recessed	Glass	White	4	0.6	113 lbs
🔒	ARS3PV	3 cu.ft.	34" x 18.5" x 19"	Freestanding	White	White	4	0.6	100 lbs
🔒	ARG3PV	3 cu.ft.	34" x 18.5" x 19"	Freestanding	Glass	White	4	0.6	105 lbs
🔒	ARS6PV	6 cu.ft.	32.5" x 23.38" x 24.38"	Freestanding	White	White	4	0.75	120 lbs
🔒	ARG6PV	6 cu.ft.	32.5" x 23.38" x 24.38"	Freestanding	Glass	White	4	0.75	125 lbs
🔒	ARS62PVBIADA	6 cu.ft.	31.75" x 23.38" x 25"	Recessed	White	White	4	0.75	125 lbs
🔒	ARG61PVBIADA	6 cu.ft.	31.75" x 23.38" x 25"	Recessed	Glass	White	4	0.75	125 lbs
🔒	ARS8PV	8 cu.ft.	50" x 23.38" x 24.38"	Freestanding	White	White	7	1.4	170 lbs
🔒	ARG8PV	8 cu.ft.	50" x 23.38" x 24.38"	Freestanding	Glass	White	7	1.4	175 lbs
🔒	ARS12PV	12 cu.ft.	61.75" x 23.38" x 24.38"	Freestanding	White	White	7	1.4	185 lbs
🔒	ARG12PV	12 cu.ft.	61.75" x 23.38" x 24.38"	Freestanding	Glass	White	7	1.4	190 lbs
🔒	ARS15PV	15 cu.ft.	72" x 23.38" x 24.38"	Freestanding	White	White	7	1.4	215 lbs
🔒	ARG15PV	15 cu.ft.	72" x 23.38" x 24.38"	Freestanding	Glass	White	7	1.4	215 lbs
🔒	ARS15PVLOCKER	15 cu.ft.	72" x 23.38" x 24.38"	Freestanding	White	White	8 lockers	1.4	250 lbs
🔒	ARG15PVLOCKER	15 cu.ft.	72" x 23.38" x 24.38"	Freestanding	Glass	White	8 lockers	1.4	250 lbs

■ All units can be combined with Accucold freezers using our stacking racks

*Specifications are subject to change without notice. For the most current information, visit accucold.com

Pharma-Lab Performance Series

Refrigerators and freezers designed and **purpose-built to meet the demanding needs for pharmacy, medication, vaccine, laboratory, and life science applications**



Advanced Temperature Control

Intelligent microprocessor digital temperature controller
Precise temperature control supports meeting the CDC/VFC storage guidelines
Improved temperature accuracy with a separate temperature probe for continuous product readout
Additional sensor reads interior air temperature for superior temperature control
Adjustable operating control range; Refrigerators: +2°C to +10°C (+35.6° F to +50° F)/ Freezers: -30°C to -15°C (-22°F to +5°F)
Digital display of the min/max temperature (updated in 15 minute increments)
Temperature viewable in Celsius or Fahrenheit to the nearest degree
Optimized forced air cooling for excellent temperature stability and uniformity with rapid recovery after door openings and product loading



Multiple Alarms For Safer Monitoring

Audible/visual alarms for high/low temperature excursions
Open door alarm will sound at different levels when the door is left ajar for times starting at 30 seconds
Condenser cleaning alarm sends audible and visual signal every six months



Fully Featured Interior For Convenient Storage

Adjustable plastic-coated wire shelves for improved circulation
Multiple shelf positions spaced at 1/2" to allow flexible adjusting to support smaller items
Interior LED lighting for an energy efficient presentation
Adjustable cycle defrost operation



Quality Design For Reliable Performance

Complete stainless steel construction inside and out for easier cleanup
Glass door models feature anti-condensate door frame heater to ensure a clear view
Self-closing door helps to protect contents, but will also remain open when the door is opened at a 90° angle for easier loading
Ships with 6" high commercial locking casters for optional installation
Factory-installed locks located towards the top of the unit for convenient security
Sealed back for a space-saving design that's easier to clean
Environmentally-friendly natural hydrocarbon refrigerant for reduced global warming potential
Meets EPA and SNAP requirements
Quality door gaskets ensure a tight seal for energy savings and improved performance

SKU: ARG23ML



6" high locking castors



Adjustable shelving



Microprocessor digital temperature control



ARS23ML



ARG49ML



AFS49ML

Technical Specifications (All Units):

Defrost type: Cycle
Thermostat: Intelligent microprocessor digital
Temperature sensor: 2 (air probe & product probe)
Adjustable operating control range: 2 to 10°C (refrigerators)/-30 to -15°C (freezers)
Shelving: White plastic coated wire
Voltage: 115V/60Hz
Cord: NEMA 5-15P
Lock: Factory-installed key lock
Refrigerant: Eco-friendly R600a (refrigerators)
 Eco-friendly R290a (Freezers)
Light: LED

Door: RHD/LHD or French (double door)
Fan: Internal
Certifications: UL, UL-C, ETL-S
Parts & Labor Warranty: 2 years
Compressor Warranty: 5 years

Included with units:
 2 keys
 User manual

Designed to operate within +/-1.5°C (Refrigerators) or +/-2.5°C (Freezers) of setpoint. Freezers with auto defrost may see temperature excursions beyond this during defrost. The published temperature variation is derived from the maximum deviation of a NTC sensor in 1 oz. vial located nearest the chamber geometric center during a 24hr test period.



Extended Warranty
 2-year parts & labor
 5-year compressor



Eco-Friendly
 Uses environmentally friendly hydrocarbons for a greener footprint



Logging Option
 Add our D2LB data logger for additional assurance



Commercially Approved
 ETL-S listed to NSF-7 commercial standards

Refrigerators

MODEL	CAPACITY	HEIGHT x WIDTH x DEPTH	DOOR(S)	CABINET	DOOR SWING	SHELVES	AMPERAGE	SHIPPING WT.
ARS23ML	23 cu.ft.	83.75" x 27.5" x 31"	(1) Stainless Steel	Stainless Steel	RHD	4	7.5	292 lbs
ARS23MLLH	23 cu.ft.	83.75" x 27.5" x 31"	(1) Stainless Steel	Stainless Steel	LHD	4	7.5	292 lbs
ARG23ML	23 cu.ft.	83.75" x 27.5" x 31"	(1) Glass	Stainless Steel	RHD	4	7.5	292 lbs
ARG23MLLH	23 cu.ft.	83.75" x 27.5" x 31"	(1) Glass	Stainless Steel	LHD	4	7.5	292 lbs
ARS49ML	49 cu.ft.	83.75" x 55.25" x 31"	(2) Stainless Steel	Stainless Steel	French	8	7.5	430 lbs
ARG49ML	49 cu.ft.	83.75" x 55.25" x 31"	(2) Glass	Stainless Steel	French	8	7.5	450 lbs

Freezers

MODEL	CAPACITY	HEIGHT x WIDTH x DEPTH	DOOR(S)	CABINET	DOOR SWING	SHELVES	AMPERAGE	SHIPPING WT.
AFS23ML	23 cu.ft.	83.75" x 27.5" x 31"	(1) Stainless Steel	Stainless Steel	RHD	4	8	292 lbs
AFS23MLLH	23 cu.ft.	83.75" x 27.5" x 31"	(1) Stainless Steel	Stainless Steel	LHD	4	8	292 lbs
AFS49ML	49 cu.ft.	83.75" x 55.25" x 31"	(2) Stainless Steel	Stainless Steel	French	8	10.5	450 lbs

*Specifications are subject to change without notice. For the most current information, visit accucold.com

Med-Lab Performance Series

Medical grade, temperature-controlled solutions to meet the demanding needs of scientific, life science, laboratory, and medical cold storage applications



Advanced Temperature Control

- Intelligent microprocessor** digital temperature controller
- Improved temperature accuracy** with a buffered probe in glycol-filled bottle to simulate temperature of stored product
- Additional sensor** reads interior air temperature for superior temperature control
- Adjustable operating control range** from +2°C to +10°C (+35.6°F to +50°F)
- Digital display** of the min/max temperature (updated in 15 minute increments)
- Temperature viewable in Celsius and Fahrenheit** to the nearest tenth of a degree
- Password protection** of control parameters beyond setpoint
- Optimized forced air cooling** for excellent temperature stability and uniformity with rapid recovery after door openings and product loading



Multiple Alarms For Safer Monitoring

- Audible/visual alarms** for high/low temperature excursions
- Open door alarm** sounds if the door is left ajar (default 4 minutes with selectable 1-99 minute range)
- Audible power failure** and sensor failure alarms
- Mute function** for alarms
- Remote alarm** contacts in rear of the unit



Fully Featured Interior For Convenient Storage

- Energy efficient LED lighting** with internal on/off switch
- Adjustable plastic-coated wire shelves** for improved circulation
- Multiple shelf positions** spaced at 1/2" to allow flexible adjusting to support smaller items
- Adjustable cycle defrost** operation
- 1/2" **probe port** included to accept user-provided monitoring equipment



Quality Design For Reliable Performance

- Hospital grade cord** with a green dot plug
- Antimicrobial sliver ion-coated handle** for added protection from germs
- Self-closing door** offers added storage protection
- User-reversible door swing** (self-closing function compatible with RHD position only)
- Environmentally-friendly** natural hydrocarbon refrigerant for reduced global warming potential
- Meets EPA and SNAP** requirements
- Convenient security** with a top-mounted door lock (2 keys included)
- Magnetic gaskets** to ensure a positive seal
- Compatible with** Pyxis®, Omnicell® and AcuDose RX®



SKU: ARG3ML



Remote Contacts



Temperature sensor in glycol bottle





New Options Available

- Keyless Locks
- Interior Lockers
- Ventilated Drawers

See Page 5



ARG6ML



ARS6MLDR



ARS8ML



ARS15MLLOCKER

Technical Specifications (All Units):

Defrost type: Cycle
Thermostat: Intelligent microprocessor digital
Temperature sensors: 2 (air probe & buffered probe)
Adjustable operating control range: +2°C to +10°C (+35.6°F to +50°F)
Temperature variation: +/- 1°C*
Shelving: Adjustable white plastic coated wire
Voltage: 115V/60 Hz
Cord: NEMA 5-15P
Lock: Factory-installed key lock
Refrigerant: Eco-friendly R600a
Light: LED
Access port: 1/2"

Door: Reversible
Fan: Internal
Certifications: ETL, ETL-C
Parts & Labor Warranty: 2 years
Compressor warranty: 5 years

Included with units:
 2 keys
 Do Not Unplug removable sticker
 Do Not Adjust Temperature removable sticker
 Do Not Store Personal Items removable sticker
 User manual

DESIGNED TO PROVIDE +/-1°C VARIATION
 The published temperature variation is derived from the maximum deviation of an NTC sensor in a 1 oz. vial located nearest the chamber geometric center during a 24 hour test period.



Extended Warranty
 2-year parts & labor
 5-year compressor



Eco-Friendly
 Uses environmentally friendly hydrocarbons for a greener footprint



Compatible
 with Pyxis®, Omnicell®, & AcuDose RX®



Logging Option
 add our D2LB data logger for additional assurance


	MODEL	CAPACITY	HEIGHT x WIDTH x DEPTH	DOOR	CABINET	SHELVES	AMPERAGE	SHIPPING WT.
🔒	ARS1ML	1 cu.ft.	21.5" x 17.5" x 19.75"	White	White	3	0.6	70 lbs
🔒	ARG1ML	1 cu.ft.	21.5" x 17.5" x 19.75"	Glass	White	3	0.6	73 lbs
🔒	ARS3ML	3 cu.ft.	34" x 18.5" x 19"	White	White	4	0.6	90 lbs
🔒	ARG3ML	3 cu.ft.	34" x 18.5" x 19"	Glass	White	4	0.6	95 lbs
🔒	ARS6ML	6 cu.ft.	32.5" x 23.38" x 24.38"	White	White	4	0.75	120 lbs
🔒	ARG6ML	6 cu.ft.	32.5" x 23.38" x 24.38"	Glass	White	4	0.75	125 lbs
🔒	ARS8ML	8 cu.ft.	50" x 23.38" x 24.38"	White	White	7	1.4	170 lbs
🔒	ARG8ML	8 cu.ft.	50" x 23.38" x 24.38"	Glass	White	7	1.4	175 lbs
🔒	ARS12ML	12 cu.ft.	61.75" x 23.38" x 24.38"	White	White	7	1.4	185 lbs
🔒	ARG12ML	12 cu.ft.	61.75" x 23.38" x 24.38"	Glass	White	7	1.4	190 lbs
🔒	ARS15ML	15 cu.ft.	72" x 23.38" x 24.38"	White	White	7	1.4	215 lbs
🔒	ARG15ML	15 cu.ft.	72" x 23.38" x 24.38"	Glass	White	7	1.4	215 lbs
🔒	ARS15MLLOCKER	15 cu.ft.	72" x 23.38" x 24.38"	White	White	8 lockers	1.4	250 lbs
🔒	ARG15MLLOCKER	15 cu.ft.	72" x 23.38" x 24.38"	Glass	White	8 lockers	1.4	250 lbs

■ All units can be combined with Accucold freezers using our stacking racks

*Specifications are subject to change without notice. For the most current information, visit accucold.com

MED2 Medical Freezers

Countertop, ADA compliant, and undercounter all-freezers designed to provide **improved temperature control** and **convenient features for pharmacy, vaccine, medical, biological, laboratory, life science, and industrial storage applications**



Advanced Temperature Control

Digital control panel for easy temperature management
Temperature viewable in Celsius or Fahrenheit to the nearest tenth of a degree
Choices for -25 or -20°C freezer storage



NIST Calibrated Monitoring

Alarm/thermometer with external readout of the current and high/low temperature
Buffered temperature probe in glycol-filled bottle
Certificate of NIST calibration included (2 years)
Audible alarm for high/low temperature excursions



User-Friendly Features

Adjustable shelves or slide-out drawers for easy storage
Hospital grade cord and plug for added assurance in high traffic settings
Front-breathing options for built-in installation under counters
3/8" access port included for additional monitoring equipment



Quality Design For Reliable Performance

Factory-installed locks front-located for convenient security
Powerful magnetic door gaskets to ensure a tight seal for energy savings and improved performance



SKU: VT65MLBIMED2ADA



Storage baskets with plastic fronts for added protection



User-friendly control panel with integrated alarm & thermostat

accucold



Technical Specifications (All Units):

Defrost type: Manual
Thermostat: Intelligent microprocessor digital
Temperature sensor: 1 (product probe)
Shelving: White plastic coated wire or slide-out plastic baskets
Voltage: 115V/60Hz
Cord: NEMA 5-15P with hospital grade plug
Lock: Factory-installed key lock
Door: RHD; factory reversible to LHD upon order
Parts & Labor Warranty: 2 years
Compressor Warranty: 5 years

Included with units:
 2 keys
 User manual
 Certificate of calibration

**DESIGNED TO PROVIDE
 +/-4°C VARIATION**
 The published temperature variation is derived from the maximum deviation of an NTC sensor in a 1 oz. vial located nearest the chamber geometric center during a 24 hour test period.



Extended Warranty
 2-year parts & labor
 5-year compressor



Logging Option
 add our D2LB data logger
 for additional assurance



NIST Certified
 Continuous temperature display
 calibrated to NIST standards in our ISO/IEC
 17025:2015 certified laboratory


SKU	CAPACITY	HEIGHT x WIDTH x DEPTH	DOOR	CABINET	SHELVES	INSTALLATION	AMPERAGE	SHIPPING WT.	CERTIFICATIONS*
FS24LMED2	1.4 cu.ft.	21.38" x 18.5" x 17.63"	White	White	1 wire	Freestanding	0.8	42 lbs	ETL, ETL-C
FS24L7MED2	1.4 cu.ft.	21.75" x 18.5" x 17.63"	White	White	1 wire	Freestanding	0.8	42 lbs	ETL, ETL-C, ETL-S
FS30LMED2	1.8 cu.ft.	24.63" x 18.5" x 17.75"	White	White	2 wire	Freestanding	0.8	53 lbs	ETL, ETL-C
FS30L7MED2	1.8 cu.ft.	26.5" x 18.5" x 17.75"	White	White	2 wire	Freestanding	0.8	53 lbs	ETL, ETL-C, ETL-S
VT65MLBIMED2	3.5 cu.ft.	33.5" x 23.63" x 23.63"	White	White	3 baskets	Built-In/Freestanding	1.3	110 lbs	ETL, ETL-C
VT65MLBIMED2ADA	3.5 cu.ft.	32" x 23.63" x 23.50"	White	White	3 baskets	Built-In/Freestanding	1.3	110 lbs	ETL, ETL-C
VT65MLBI7MED2	3.5 cu.ft.	33.75" x 23.63" x 23.5"	White	White	3 baskets	Built-In/Freestanding	1.3	110 lbs	ETL, ETL-C, ETL-S
VT65MLBI7MED2ADA	3.5 cu.ft.	32" x 23.63" x 23.50"	White	White	3 baskets	Built-In/Freestanding	1.3	110 lbs	ETL, ETL-C, ETL-S

All units can be combined with Accucold refrigerators using our stacking racks

*Specifications are subject to change without notice. For the most current information, visit accucold.com

ACR Pharmacy Refrigerators

ADA compliant and full-sized **all-refrigerators** designed for the **reliable storage of pharmaceuticals, vaccines, and medical/biological substances**



Features

- Commercially listed** to NSF-7 standards and UL-471
- LED lighting** for a bright but efficient interior view
- Double pane tempered glass door** to protect stored items (select units)
- Stainless steel interior construction** on full-sized units for added durability and easier cleanup
- Self-closing door** for added protection (select units)
- Factory-installed lock** to provide added security
- Fan-cooled interior** ensures even temperature distribution
- Digital control panel** with a +2 to 8°C temperature range
- Hospital grade cord** and plug for added assurance in high traffic settings



Drawer option available for easy installation inside the unit.



SKU: ACR45L




SKU: ACR1718RH

SKU	CAPACITY	HEIGHT x WIDTH x DEPTH	DOOR	CABINET	DOOR SWING	SHELVES	AMPERAGE	SHIPPING WT.	CERTIFICATIONS
🔒 ACR45L	2.6 cu.ft.	31.5" x 19.5" x 22.75"	White	White	Reversible	3 wire	1	89 lbs	ETL, ETL-C, ETL-S
🔒 ACR1415RH	14 cu.ft.	66" x 23.63" x 24.25"	Glass	Stainless Steel	RHD	4 molded plastic	4.4	177 lbs	ETL, ETL-C, ETL-S
🔒 ACR1415LH	14 cu.ft.	66" x 23.63" x 24.25"	Glass	Stainless Steel	LHD	4 molded plastic	4.4	177 lbs	ETL, ETL-C, ETL-S
🔒 ACR1718RH	17 cu.ft.	79.25" x 27.63" x 24.75"	Glass	Stainless Steel	RHD	6 molded plastic	5	220 lbs	ETL, ETL-C, ETL-S
🔒 ACR1718LH	17 cu.ft.	79.25" x 27.63" x 24.75"	Glass	Stainless Steel	LHD	6 molded plastic	5	220 lbs	ETL, ETL-C, ETL-S

Stacked Refrigerator-Freezer Combinations

Pharma-Vac refrigerators stacked with MED2 freezers create **all-purpose combinations** that provide reliable storage in slimmer fitting footprints



Features

- Easy assembly** with our stacking rack
- Takes up less than 24" of floor space** for an easier fit in facilities
- Fully featured** with our PV and MED2 packages
- Access ports included** with both units
- External temperature readout** as recommended by the CDC
- Independent refrigerator and freezer operation** for optimum performance
- High/low temperature alarms** installed in both units
- Locking doors** with separate keys
- Extended warranty** with 2 years parts & labor/5 years compressor



StackRack lets you combine virtually any two units



SKU: ARS6PV-VT65MLSTACKMED2



SKU: ARG8PV-FS24LSTACKMED2

SKU	CAPACITY (REF/FRZ)	HEIGHT x WIDTH x DEPTH	DOOR (REF/FRZ)	CABINET	SHELVES (REF/FRZ)	AMPERAGE (REF/FRZ)	CERTIFICATIONS
ARS3PV-ADA305AFSTACK	3/2.47 cu.ft.	66" x 19.75" x 22.5"	White	White	4/3	.6/1.1	ETL, ETL-C
ARG3PV-ADA305AFSTACK	3/2.47 cu.ft.	66" x 19.75" x 22.5"	Glass/White	White	4/3	.6/1.1	ETL, ETL-C
ARS6PV-VT65MLSTACKMED2	6/3.2 cu.ft.	68" x 23.63" x 24.38"	White	White	4/3	.75/1.3	ETL, ETL-C
ARG6PV-VT65MLSTACKMED2	6/3.2 cu.ft.	68" x 23.63" x 24.38"	Glass/White	White	4/3	.75/1.3	ETL, ETL-C
ARS8PV-FS24LSTACKMED2	8/1.4 cu.ft.	73.5" x 23.38" x 24.38"	White	White	7/1	1.4/.8	ETL, ETL-C
ARG8PV-FS24LSTACKMED2	8/1.4 cu.ft.	73.5" x 23.38" x 24.38"	Glass/White	White	7/1	1.4/.8	ETL, ETL-C
ARS8PV-FS30LSTACKMED2	8/1.8 cu.ft.	76.5" x 23.38" x 24.38"	White	White	7/2	1.4/.8	ETL, ETL-C
ARG8PV-FS30LSTACKMED2	8/1.8 cu.ft.	76.5" x 23.38" x 24.38"	Glass/White	White	7/2	1.4/.8	ETL, ETL-C

* All units can be manufactured with a left hand door swing. Add option "LHD" upon order
 ■ Additional combinations available by special order

*Specifications are subject to change without notice. For the most current information, visit accucold.com




PATIENT CARE/ HOSPITAL/BREAK ROOM

From basic refrigeration and towel warming cabinets to breast milk storage refrigeration for hospitals or corporate offices, Accucold's additional product lines help to serve a variety of facility needs.

PURETHERM Blanket Warmers

Designed and constructed in the USA to keep towels, blankets, and other linens properly warmed in hospitals, clinics, spas, and other facilities



Features

- Passive heating system** (PTHC series) to minimize lint build up for safe, effective operation
- Stainless steel or hammered aluminum interior** for easy cleaning
- Double pane tempered glass door** provides a full view of blanket inventory
- Digital thermostats** ensure intuitive temperature management
- Large LED readout** of the temperature
- Operates between 90 and 140°F**
- Factory-installed lock** with two keys included
- Silent operation** (PTHC series) with no moving parts
- Automatic turnoff** prevents overheating
- Flush design** with a sealed back
- Optional dollies** available for easier mobility
- Fully assembled** in the USA



Stainless steel construction for lasting durability



SKU: PTHC65G




SKU: PTHC155G

	SKU	CAPACITY	HEIGHT x WIDTH x DEPTH	DOOR	CABINET	DOOR SWING	AMPERAGE	SHIPPING WT.	CERTIFICATIONS
🔒	PHC51G	2 cu.ft.	29.5" x 23.38" x 19.88"	Glass	Gray	RHD	5.5	145 lbs.	ETL, ETL-C
🔒	PTHC65G	4.57 cu.ft.	36" x 23.63" x 23.5"	Glass	Stainless Steel	RHD	5.5	145 lbs.	ETL, ETL-C
🔒	PTHC65GLHD	4.57 cu.ft.	36" x 23.63" x 23.5"	Glass	Stainless Steel	LHD	5.5	145 lbs.	ETL, ETL-C
🔒	PTHC125G	10.57 cu.ft.	60" x 23.63" x 23.75"	Glass	Stainless Steel	RHD	5.5	175 lbs.	ETL, ETL-C
🔒	PTHC155G	12.17 cu.ft.	71.5" x 23.63" x 23.75"	Glass	Black	RHD	5.5	200 lbs.	ETL, ETL-C
🔒	PTHC155GLHD	12.17 cu.ft.	71.5" x 23.63" x 23.75"	Glass	Black	LHD	5.5	200 lbs.	ETL, ETL-C
🔒	PTHC155GCSS	12.17 cu.ft.	71.5" x 23.63" x 23.75"	Glass	Stainless Steel	RHD	5.5	215 lbs.	ETL, ETL-C
🔒	PTHC155GCSSLHD	12.17 cu.ft.	71.5" x 23.63" x 23.75"	Glass	Stainless Steel	LHD	5.5	215 lbs.	ETL, ETL-C

*Specifications are subject to change without notice. For the most current information, visit accucold.com

JCAHO Compliant Nutrition Centers

Commercially approved refrigerators and freezers specifically designed for the storage of sealed food and beverages in a variety of facilities including break rooms, waiting areas, nurses' stations, and more



Features

- Factory-installed NIST calibrated thermometer with continuous display
- 2-year certificate of calibration included for your records
- Commercially approved to NSF-7 standards
- Audible high/low temperature alarms
- Factory-installed keyed locks
- Adjustable shelving for easier storage
- Matching refrigerator and freezer sets available
- Range of sizes including countertop, undercounter, ADA compliant, and full capacity



SKU: SCR215LNZ



SKU: FF31L7NZ



SKU: SCR600BGLBINZ



SKU: FFAR121SSNZ

Refrigerators

SKU	CAPACITY	HEIGHT x WIDTH x DEPTH	DOOR	CABINET	SHELVES	INSTALLATION	AMPERAGE	SHIPPING WT.	CERTIFICATIONS
FFAR25L7NZ	1.7 cu.ft.	22.75" x 17" x 19.13"	White	White	2	Freestanding	1.3	60 lbs	ETL, ETL-C, ETL-S
SCR215LNZ	1.7 cu.ft.	22.75" x 17" x 18.75"	Glass	White	2	Freestanding	1.3	60 lbs	ETL, ETL-C, ETL-S
SCR312LNZ	2.5 cu.ft.	28" x 17" x 19.25"	Glass	Black	4	Freestanding	1.3	100 lbs	ETL, ETL-C, ETL-S
SCR486LNZ	3.35 cu.ft.	35.5" x 19" x 19.5"	Glass	Black	4	Freestanding	1.3	100 lbs	ETL, ETL-C, ETL-S
FF6LWBI7NZ	5.5 cu.ft.	33.5" x 23.63" x 23.5"	White	White	3	Built-in/Freestanding	1.3	110 lbs	ETL, ETL-C, ETL-S
FF6LWBI7NZADA	5.5 cu.ft.	32" x 23.63" x 23.5"	White	White	3	Built-in/Freestanding	1.3	110 lbs	ETL, ETL-C, ETL-S
SCR600GLBINZ	5.5 cu.ft.	33.5" x 23.63" x 23.5"	Glass	White	4	Built-in/Freestanding	1.8	122 lbs	UL, UL-C, ETL-S
SCR600GLBINZADA	5.5 cu.ft.	32" x 23.63" x 23.5"	Glass	White	4	Built-in/Freestanding	1.8	122 lbs	UL, UL-C, ETL-S
FFAR12WNZ	10.1 cu.ft.	60.88" x 23.63" x 22.75"	White	White	4	Freestanding	1.7	155 lbs	ETL, ETL-C, ETL-S
FFAR121SSNZ	10.1 cu.ft.	60.88" x 23.63" x 22.75"	Stainless Steel	Stainless Steel	4	Freestanding	1.7	155 lbs	ETL, ETL-C, ETL-S

Freezers

SKU	CAPACITY	HEIGHT x WIDTH x DEPTH	DOOR	CABINET	SHELVES	INSTALLATION	AMPERAGE	SHIPPING WT.	CERTIFICATIONS
FS24L7NZ	1.4 cu.ft.	22.25" x 18.5" x 17.63"	White	White	1	Freestanding	0.8	55 lbs	ETL, ETL-C, ETL-S
SCFU386NZ	2 cu.ft.	26.25" x 23.38" x 20.63"	Glass	Gray	2	Freestanding	4	95 lbs	ETL, ETL-C, ETL-S

*Specifications are subject to change without notice. For the most current information, visit accucold.com

Break Room Refrigerator-Freezers

Variety of choices equipped with NIST calibrated thermometers that provide a continuous digital display of the current and high/low temperature



Features

- Factory-installed NIST calibrated thermometer**
- 2-year certificate of calibration** included for your records
- Audible high/low temperature alarm**
- Digital temperature readout** for each section (refrigerator only for undercounter choices)
- Frost-free operation** in full-sized units
- Adjustable** shelving and door storage
- Interior lighting** in the fresh food compartment
- User-reversible door swing** to best accommodate any setting



SKU: BKRF663B



SKU: BKRF661



SKU: BKRF18SSCP with antimicrobial copper handles



SKU: BKRF21W

SKU	CAPACITY	HEIGHT x WIDTH x DEPTH	DOOR	CABINET	SHELVES	INSTALLATION	AMPERAGE	SHIPPING WT.	CERTIFICATIONS
BKRF661	5.1 cu.ft	32.63" x 23.63" x 23"	White	White	3 glass	Freestanding	1.3	110 lbs	UL, UL-C
BKRF661BIADA	5.1 cu.ft	32" x 23.63" x 23"	White	White	3 glass	Built-in/Freestanding	1.3	110 lbs	UL, UL-C
BKRF663B	5.1 cu.ft	32.63" x 23.63" x 23"	Black	Black	3 glass	Freestanding	1.3	110 lbs	UL, UL-C
BKRF663BBIADA	5.1 cu.ft	32" x 23.63" x 23"	Black	Black	3 glass	Built-in/Freestanding	1.3	110 lbs	UL, UL-C
BKRF1118W	10.3 cu.ft	60.63" x 23.63" x 26"	White	White	3	Freestanding	1.6	105 lbs	UL, UL-C
BKRF1119BL	10.3 cu.ft	60.63" x 23.63" x 26"	Black	Black	3	Freestanding	1.6	105 lbs	UL, UL-C
BKRF1159SS	10.3 cu.ft	60.63" x 23.63" x 26"	Black	Stainless Steel	3	Freestanding	1.6	170 lbs	UL, UL-C
BKRF15W	14.75 cu.ft	61.63" x 28" x 29.38"	White	White	3	Freestanding	4.5	182 lbs	CSA
BKRF18W	18 cu.ft	67.63" x 29.63" x 30"	White	White	3	Freestanding	6	220 lbs	CSA
BKRF18B	18 cu.ft	67.63" x 29.63" x 30"	Black	Black	3	Freestanding	6	220 lbs	CSA
BKRF18SS	18 cu.ft	67.63" x 29.63" x 30"	Black	Stainless Steel	3	Freestanding	6	220 lbs	CSA
BKRF21W	20.4 cu.ft	70.5" x 29.63" x 31.75"	White	White	4	Freestanding	6	230 lbs	CSA
BKRF21B	20.4 cu.ft	70.5" x 29.63" x 31.75"	Black	Black	4	Freestanding	6	230 lbs	CSA
BKRF21SS	20.4 cu.ft	70.5" x 29.63" x 31.75"	Black	Stainless Steel	4	Freestanding	6	230 lbs	CSA

*Most units can be ordered with combination or key locks

*Specifications are subject to change without notice. For the most current information, visit accucold.com

MOMCUBE Series Breast Milk Refrigeration

All-refrigerators specially designed for the unique needs of breast milk, with interior locking compartments to allow nursing mothers private storage



Countertop



ADA



Undercounter



Full Sized



Purpose-Built for Breast Milk Storage

Designed to operate between 33 and 40°F for the proper storage of human milk
Individual locking compartments for secure, separate storage for ("LK" Series)
Combination locks allow each mother her own access system
Single key override included for staff to manage as needed
Audible/visual alarms for high/low temperature excursions
Open door alarm offers added assurance in office or hospital settings



Advanced Temperature Performance

Intelligent microprocessor digital temperature controller
Improved accuracy with a buffered probe in glycol-filled bottle to simulate temperature of stored bottles
Additional sensor reads interior air temperature for superior temperature control
Digital display of the min/max temperature (updated in 15 minute increments)
Password protection of control parameters beyond setpoint
Optimized forced air cooling for excellent temperature stability



Fully Featured Interior For Convenient Storage

Energy efficient LED lighting with internal on/off switch
1/2" probe port included to accept user-provided monitoring equipment



Quality Design For Reliable Performance

Hospital grade cord with a green dot plug
Antimicrobial silver ion-coated handle for added protection from germs
Self-closing door offers added storage protection
Environmentally friendly natural hydrocarbon refrigerant for reduced global warming potential
Meets EPA and SNAP requirements
Convenient security with a top-mounted front door lock (2 keys included)
Powerful magnetic door gaskets to ensure a tight seal for energy savings and improved performance



SKU: ARS8MLMCLK

MOMCUBE Freezers:

Provide additional solutions for longer term storage. These units come with a more basic featured package.



Interior locking compartments allow private storage



Designed to maintain recommended temperature for breast milk storage



Technical Specifications (Refrigerators*):

Defrost type: Cycle
Thermostat: Intelligent microprocessor digital
Temperature sensors: 2 (air probe & buffered probe)
Adjustable operating control range: 33 to 40°F
Temperature variation: +/-1°C*
Shelving: Interior locking compartments or adjustable white plastic coated wire
Voltage: 115V/60Hz
Cord: NEMA 5-15P
Lock: Factory-installed key lock + individual combination locks with single key override ("LK" Series)
Refrigerant: Eco-friendly R600a
Light: LED

Access Port: 1/2"
Door: Reversible
Fan: Internal
Certifications: ETL, ETL-C
Parts & Labor Warranty: 2 years
Compressor Warranty: 5 years

Included with units:
 2 keys for outer door
 2 keys to override combination locks
 User manual

DESIGNED TO PROVIDE +/-1°C VARIATION
 The published temperature variation is derived from the maximum deviation of an NTC sensor in a 1 oz. vial located nearest the chamber geometric center during a 24 hour test period.



Extended Warranty
 2-year parts & labor
 5-year compressor



Eco-Friendly
 Uses environmentally friendly hydrocarbons for a greener footprint



Logging Option
 Add our D2LB data logger for additional assurance



LOCKERS
 Individual storage compartments with combination locks

Refrigerators

MODEL	CAPACITY	HEIGHT x WIDTH x DEPTH	INSTALLATION	DOOR	CABINET	INTERIOR	AMPERAGE	SHIPPING WT.
ARS1MLMC	1 cu.ft.	21.5" x 17.5" x 19.75"	Freestanding	White	White	2 shelves	0.6	70 lbs
ARS3MLMC	3 cu.ft.	34" x 18.5" x 19"	Freestanding	White	White	4 shelves	0.6	90 lbs
ARS32MLMCBIADA	2.83 cu.ft.	32" x 23.38" x 22.63"	Recessed	White	White	4 shelves	0.6	113 lbs
ARS6MLMC	6 cu.ft.	32.5" x 23.38" x 24.63"	Freestanding	White	White	4 shelves	0.75	120 lbs
ARS6MLMCLK	6 cu.ft.	32.5" x 23.38" x 24.63"	Freestanding	White	White	2 locking compartments	0.75	150 lbs
ARS62MLMCBIADA	6 cu.ft.	32" x 23.38" x 25"	Recessed	White	White	4 shelves	0.75	150 lbs
ARS62MLMCBIADADLK	6 cu.ft.	32" x 23.38" x 25"	Recessed	White	White	2 locking compartments	0.75	150 lbs
ARS8MLMC	8 cu.ft.	50" x 23.38" x 24.38"	Freestanding	White	White	7 shelves	1.4	170 lbs
ARS8MLMCLK	8 cu.ft.	50" x 23.38" x 24.38"	Freestanding	White	White	2 locking compartments	1.4	210 lbs
ARS12MLMC	12 cu.ft.	61.75" x 23.38" x 24.38"	Freestanding	White	White	7 shelves	1.4	185 lbs
ARS12MLMCLK	12 cu.ft.	61.75" x 23.38" x 24.38"	Freestanding	White	White	6 locking compartments	1.4	230 lbs
ARS15MLMC	15 cu.ft.	72" x 23.38" x 24.38"	Freestanding	White	White	7 shelves	1.4	215 lbs
ARS15MLMCLK	15 cu.ft.	72" x 23.38" x 24.38"	Freestanding	White	White	8 locking compartments	1.4	300 lbs

Freezers

MODEL	CAPACITY	HEIGHT x WIDTH x DEPTH	INSTALLATION	DOOR	CABINET	INTERIOR	AMPERAGE	SHIPPING WT.
FS24LMC	1.4 cu.ft.	20" x 18.5" x 17.63"	Freestanding	White	White	1 wire shelf	0.8	37 lbs
FS30LMC	1.8 cu.ft.	24.5" x 18.5" x 17.75"	Freestanding	White	White	2 wire shelves	0.8	42 lbs
ALFZ36MC	2.68 cu.ft.	32" x 19.75" x 22.5"	Recessed	White	White	3 wire shelves	0.6	70 lbs
ALFZ36LMCTBC	2.68 cu.ft.	32" x 19.75" x 22.5"	Recessed	White	White	3 wire shelves	0.6	70 lbs
VT65MLBIMCADA	3.5 cu.ft.	32" x 23.63" x 23.5"	Recessed	White	White	3 baskets	1.3	110 lbs

* For complete freezer technical specifications, see our website

Compact & Mid-Sized Microwaves

Commercially approved units and basic countertop microwaves for a variety of heating applications



Features

- Black, white, or stainless steel** finish options
- Commercially approved** choices listed to NSF-4 standards
- 115V operation** to work with standard outlets
- Digital control** pads for easy use
- Savable settings** for faster one-touch cooking
- Can be ordered pre-mounted** to virtually any Accucold or Summit refrigerator



Intuitive digital controls allow you to store settings for faster cooking




SKU: SCM1000SS

SKU	CAPACITY	HEIGHT x WIDTH x DEPTH	FINISH	AMPERAGE	SHIPPING WT.	WATTAGE	CERTIFICATIONS
SCM853	.7 cu.ft.	11" x 18.88" x 14.25"	Stainless Steel	12	32 lbs	900W	UL, UL-C
SM901WH	.7 cu.ft.	10.25" x 17.75" x 13"	White	15	37 lbs	700W	UL, UL-C
SM902BL	.7 cu.ft.	10.25" x 17.75" x 13"	Black	15	37 lbs	700W	UL, UL-C
SM1102WH	1 cu.ft.	11.75" x 21.13" x 15.5"	White	20	37 lbs	1000W	UL, UL-C
SCM1000SS	.9 cu.ft.	12.25" x 20.25" x 14.75"	Stainless Steel	20	41 lbs	1000W	UL, UL-C, ETL-S

Icemakers & Ice/Water Dispensers

Available in a range of sizes, including countertop, undercounter, and ADA compliant for use under lower counters



Features

- Tabletop ice/water dispenser ideal for hospital waiting areas
- Commercially approved choices listed to NSF-12 standards
- Drain-free choice (BIM25) for versatile installation
- Choices available in a lower 32" height for use under ADA compliant counters
- ENERGY STAR certified choices included
- High capacity production on select models
- Available with an internal pump or gravity drain system



SKU: AIWD282

SKU: BIM100

SKU: BIM25

Ice & Water Dispenser

SKU	CAPACITY	ICE TYPE	HEIGHT x WIDTH x DEPTH	DOOR	CABINET	AMPERAGE	SHIPPING WT.	CERTIFICATIONS
AIWD282	282 lbs./day	Clear nugget	28" x 16.5" x 24.88"	Stainless Steel	Stainless Look	3.8	177 lbs	ETL, ETL-C, UL-C

Clear Icemakers (Drain Required)

SKU	CAPACITY	ICE TYPE	HEIGHT x WIDTH x DEPTH	DOOR	CABINET	AMPERAGE	SHIPPING WT.	CERTIFICATIONS
BIM44G	50 lbs./day	Clear cube	33" x 14.5" x 23.5"	Stainless Steel	Black	4.8	105 lbs	UL, UL-C, NSF
BIM44GADA	50 lbs./day	Clear cube	32.38" x 14.5" x 23.5"	Stainless Steel	Black	4.8	105 lbs	UL, UL-C, NSF
BIM100	100 lbs./day	Clear cube	34.5" x 20.5" x 22.25"	Black	Stainless Steel	3.8	120 lbs	ETL, ETL-C, ETL-S
BIM100ADA	100 lbs./day	Clear cube	31.5" x 20.5" x 22.25"	Black	Stainless Steel	3.8	120 lbs	ETL, ETL-C, ETL-S

Crescent Icemakers (No Drain Required)

SKU	CAPACITY	ICE TYPE	HEIGHT x WIDTH x DEPTH	DOOR	CABINET	AMPERAGE	SHIPPING WT.	CERTIFICATIONS
BIM25	12 lbs./day	Cloudy crescent	24.75" x 15" x 18.25"	Stainless Steel	Black	1.7	53 lbs	UL, UL-C

*Specifications are subject to change without notice. For the most current information, visit accucold.com

Portable Refrigerator/Freezers

Ensure proper cold chain management with mobile refrigeration designed for easy travel and reliable storage



Features

- Includes a **12V power cord** to connect to vehicle sockets
- 110V AC/DC power adapter** included for alternate use
- Factory-installed casters** for smooth mobility (most models)
- Flexible strap handle or trolley** allows easy steering
- Digital control panel** ensures easy operation
- Compressor-based cooling system** for more stability
- Can operate as a refrigerator or freezer** (temperature range: 2-8°C/-12°C or -15°C)



SKU: SPRF86M2



SKU: SPRF26M

SKU	CAPACITY	HEIGHT x WIDTH x DEPTH	FINISH	PORTABLE FEATURES	POWER	SHIPPING WEIGHT (LBS.)
SPRF26	.88 cu.ft.	17" x 21.25" x 14.25"	Gray	Handles	12/24V DC/115V	30
🔒 SPRF26M	.88 cu.ft.	17" x 21.25" x 14.25"	Gray	Collapsible trolley	12/24V DC/115V	32
🔒 SPRF36M2	1.06 cu.ft.	24.25" x 21.63" x 14.75"	Gray	Factory-installed casters w/handle	12/24V DC/115V	50
🔒 SPRF86M2	2.8 cu.ft.	25.88" x 22" x 17.5"	Gray	Factory-installed casters w/handle	12/24V DC/115V	70




General Purpose & NSF-7 Refrigeration

Broad selection of compact, undercounter, ADA compliant, and full-sized refrigerators and freezers for a wide variety of basic storage, including choices listed to NSF-7 commercial standards.

Healthcare & General Purpose Refrigerators & Freezers

Wide selection of general purpose options ideal for non-critical storage in a variety of healthcare applications including nursing stations, waiting rooms, recovery areas, and more



Features

- Factory-installed front locks** for keyed security
- 32" high choices** for use under lower ADA compliant counters
- Commercially listed choices** to meet NSF-7 and JCAHO guidelines
- Glass and solid door models** from compact to full-sized
- Temperature readouts** included on many models
- Matching** refrigerator and freezer sets
- Customizable** to meet your exact specifications



All-Refrigerators

SKU	CAPACITY	HEIGHT x WIDTH x DEPTH	DOOR	CABINET	SHELVES	INSTALLATION	AMPERAGE	SHIPPING WT.	CERTIFICATIONS
FFAR23LGP	1.7 cu.ft.	20.75" x 18.75" x 17.5"	White	White	1 Wire	Countertop	1	40	ETL, ETL-C
FF28LWHGP	2.4 cu.ft.	24.88" x 18.5" x 17.75"	White	White	1 Wire	Countertop	1.3	38	ETL, ETL-C
SCR314LDTGP	2.1 cu.ft.	26" x 18.88" x 17"	Glass	Black	2 Wire	Countertop	1.3	56	ETL, ETL-C, ETL-S
ADA404REF	3.53 cu.ft.	32" x 19.75" x 22.5"	White	White	3 wire	Built-in/Freestanding	1.1	70	ETL, ETL-C, ETL-S
FF511L	4.1 cu.ft.	34" X 19.25" X 21.38"	White	White	3 wire	Freestanding	1.4	60	ETL, ETL-C
SCR600GLBIADAGP	5.5 cu.ft.	32.38" x 23.63" x 23.5"	Glass	White	4 Glass	Built-in/Freestanding	1.8	105	UL, UL-C, ETL-S
FF7LWBIADAGP	5.5 cu.ft.	32.25" X 23.63" X 23.5"	White	White	3 Wire	Built-in/Freestanding	1.4	105	UL, UL-C, ETL-S
FFAR10GP	10 cu.ft.	55.63" x 23.63" x 23"	White	White	4 Wire	Freestanding	1.5	140	UL, UL-C

All Freezers


SKU	CAPACITY	HEIGHT x WIDTH x DEPTH	DOOR	CABINET	SHELVES	INSTALLATION	AMPERAGE	SHIPPING WT.	CERTIFICATIONS
FS24LGP	1.4 cu.ft.	20" x 18.5" x 17.63"	White	White	1 Wire	Countertop	0.8	40	ETL, ETL-C
FS30LGP	1.8 cu.ft.	24.63" x 18.5" x 17.75"	White	White	2 Wire	Countertop	0.8	45	ETL, ETL-C
ADA305AF	2.47 cu.ft.	32" x 19.75" x 22.5"	White	White	3 pull-out baskets	Built-in/Freestanding	0.6	70	ETL, ETL-C, ETL-S
VT65MLGP	3.5 cu.ft.	33.5" x 23.63" x 23.5"	White	White	3 Drawers	Freestanding	1.3	110	UL
VT65MLBIADAGP	3.5 cu.ft.	32.25" x 23.63" x 23.5"	White	White	4 Drawers	Freestanding	1.3	110	UL, UL-C

Refrigerator-Freezers

SKU	CAPACITY	HEIGHT x WIDTH x DEPTH	DOOR	CABINET	SHELVES	INSTALLATION	AMPERAGE	SHIPPING WT.	CERTIFICATIONS
ADA302RFZ	2.68 cu.ft.	32" x 19.75" x 22.5"	White	White	3 Glass	Built-in/Freestanding	1.1	70	ETL, ETL-C

PLUS2 Series Refrigerators & Freezers

All-refrigerators and matching all-freezers equipped with NIST calibrated thermometers to provide an external display in general purpose storage applications



Advanced Temperature Control

Digital control panel for easy temperature management
Temperature viewable in Celsius or Fahrenheit to the nearest tenth of a degree



NIST Calibrated Monitoring

Alarm/thermometer with continuous external readout of the current and high/low temperature
Buffered temperature probe in glycol-filled bottle
Certificate of NIST calibration included (2-years)
Audible alarm for high/low temperature excursions



User-Friendly Features

Adjustable shelves or slide-out drawers for easy storage
Hospital grade cord and plug for added assurance in high traffic settings
Front-breathing options for built-in installation under counters
Matching refrigerator and freezer sets
Variety of sizes including countertop, undercounter, ADA compliant, and larger capacity



Quality Design For Reliable Performance

Factory-installed locks front-located for convenient security
Fan-forced cooling to promote even temperature distribution
Powerful magnetic door gaskets to ensure a tight seal for energy savings and improved performance



SKU: FF6LWBIPLUS2ADA

Refrigerators

SKU	CAPACITY	HEIGHT x WIDTH x DEPTH	FINISH	SHELVES	INSTALLATION	AMPERAGE	SHIPPING WT.	CERTIFICATIONS
FF28LWHPLUS2	2.4 cu.ft.	24.5" x 18.5" x 17.75"	White	1	Freestanding	1.3	38 lbs.	ETL, ETL-C
FF6LWPLUS2	5.5 cu.ft.	33.5" x 23.63" x 23.5"	White	3	Freestanding	1.3	105 lbs.	ETL, ETL-C
FF6LWBIPLUS	5.5 cu.ft.	33.5" x 23.63" x 23.5"	White	3	Built-in/Freestanding	1.3	105 lbs.	ETL
FF6LWBIPLUS2ADA	5.5 cu.ft.	32" x 23.63" x 23.5"	White	3	Built-in/Freestanding	1.3	105 lbs.	ETL-C
FFAR10PLUS2	10.1 cu.ft.	56" x 23.63" x 23"	White	4	Freestanding	1.5	145 lbs.	ETL, ETL-C

Freezers

SKU	CAPACITY	HEIGHT x WIDTH x DEPTH	FINISH	SHELVES	INSTALLATION	AMPERAGE	SHIPPING WT.	CERTIFICATIONS
FS24LPLUS2	1.4 cu.ft.	19.63" x 18.5" x 17.63"	White	1	Freestanding	0.8	37 lbs.	ETL, ETL-C
FS30LPLUS2	1.8 cu.ft.	24.5" x 18.5" x 17.75"	White	2	Freestanding	0.8	53 lbs.	ETL, ETL-C


Refrigerator-Freezers

SKU	CAPACITY	HEIGHT x WIDTH x DEPTH	FINISH	SHELVES	INSTALLATION	AMPERAGE	SHIPPING WT.	CERTIFICATIONS
CP351WLLF2PLUS2	2.9 cu.ft.	32.88" x 18.5" x 20.25"	White	1	Freestanding	1.4	93 lbs.	ETL, ETL-C
CP351WLLF2PLUS2ADA	2.9 cu.ft.	32" x 18.5" x 20.25"	White	1	Freestanding	1.4	93 lbs.	ETL, ETL-C

More choices available. See our website for our complete selection

Low Temperature Freezers

From laboratory research to veterinary storage, Accucold's low temperature freezers serve a variety of storage functions in a range of sizes to fit any facility's space



Low Temperature Storage

Digital control panel for easy temperature management
Temperature viewable in Celsius or Fahrenheit to the nearest tenth of a degree
Choices for -45, -35, or -30°C operation
NIST calibrated thermometers included on VT/VLT series



Variety of Sizes & Styles

Front-opening choices available in standard undercounter & ADA compliant height
Chest style units from 28" to 86" wide



User-Friendly Features

Counterbalanced hinges to ensure safer loading
Easy-grip handle (most models)
Stainless steel corner protectors (VT/VLT series) for added durability
Casters included on EL series



Quality Design For Reliable Performance

Factory-installed locks for convenient security
Extra thick insulation for improved efficiency
Powerful magnetic door gaskets to ensure a tight seal for energy savings and improved performance
Commercially approved to NSF-7 standards




SKU: EL11LT

	SKU	CAPACITY	HEIGHT x WIDTH x DEPTH	FINISH	STYLE	SET POINT RANGE	AMPERAGE	SHIPPING WT.	CERTIFICATIONS
🔒	EL11LT	4.8 cu.ft.	33.5" x 28.5" x 27.25"	White	Chest	-45 to 0°C	7.2	115 lbs.	UL, UL-C, UL-S
🔒	EL21LT	8.4 cu.ft.	33.5" x 41.5" x 27.25"	White	Chest	-45 to 0°C	7.2	135 lbs.	UL, UL-C, UL-S
🔒	EL31LT	11.1 cu.ft.	33.5" x 51.38" x 27.25"	White	Chest	-45 to 0°C	7.2	157 lbs.	UL, UL-C, UL-S
🔒	EL41LT	12.8 cu.ft.	33.5" x 59.25" x 27.25"	White	Chest	-45 to 0°C	7.2	177 lbs.	UL, UL-C, UL-S
🔒	EL51LT	15.5 cu.ft.	33.5" x 67" x 27.25"	White	Chest	-45 to 0°C	7.1	220 lbs.	UL, UL-C, UL-S
🔒	VLT650	3.5 cu.ft.	33.5" x 23.63" x 23.5"	White	Front-Opening	-35 to -15°C	1.5	135 lbs.	UL, UL-C, UL-S
🔒	VLT650ADA	3.5 cu.ft.	32" x 23.63" x 23.5"	White	Front-Opening	-35 to -15°C	1.5	135 lbs.	UL, UL-C, UL-S
🔒	VLT850	9 cu.ft.	33" x 45.75" x 24"	White	Front-Opening	-35 to -15°C	1.4	110 lbs.	UL, UL-C, UL-S
🔒	VLT1250	14.4 cu.ft.	33.25" x 55.5" x 28.63"	White	Chest	-35 to -15°C	1.49	235 lbs.	UL, UL-C, UL-S
🔒	VLT1750	18 cu.ft.	33.25" x 66.75" x 28.63"	White	Chest	-35 to -15°C	1.4	175 lbs.	UL, UL-C, UL-S
🔒	VLT2250	24.8 cu.ft.	32" x 85.25" x 31.5"	White	Chest	-35 to -15°C	3.1	230 lbs.	UL, UL-C, UL-S
🔒	VT65	3.5 cu.ft.	33.5" x 23.63" x 23.5"	White	Front-Opening	-30 to -15°C	1.5	105 lbs.	UL, UL-C, UL-S
🔒	VT65ADA	3.5 cu.ft.	32" x 23.63" x 23.5"	White	Front-Opening	-30 to -15°C	1.5	105 lbs.	UL, UL-C, UL-S
🔒	VT85	9 cu.ft.	33" x 45.75" x 24"	White	Chest	-30 to -15°C	1.4	110 lbs.	UL, UL-C, UL-S
🔒	VT125	14.1 cu.ft.	33.25" x 55.5" x 28.63"	White	Chest	-30 to -15°C	1.49	250 lbs.	UL, UL-C, UL-S
🔒	VT175	18 cu.ft.	33.25" x 66.75" x 28.63"	White	Chest	-30 to -15°C	1.4	175 lbs.	UL, UL-C, UL-S
🔒	VT225	24.8 cu.ft.	32" x 85.25" x 31.5"	White	Chest	-30 to -15°C	3.1	240 lbs.	UL, UL-C, UL-S

*Specifications are subject to change without notice. For the most current information, visit accucold.com

Commercial Chest Freezers

Our **commercially listed chest series** includes **frost-free refrigerators and freezers**, as well as manual defrost chest freezers to meet a variety of **general purpose storage needs**



Features

Stainless steel corner protectors for added durability in high traffic areas

Commercially approved to NSF-7 standards

Counterbalanced hinges to ensure safer loading

Easy-grip handle (most models)

Factory-installed keyed locks (most models)

Choice of frost-free operation or static manual defrost

Digital controls with external temperature readout on frost-free EQ series

Pre-installed probe hole included on frost-free EQ series

Front-located dial controls on manual defrost SCFM series



Frost-free operation adds convenience to a design previously limited to manual defrost



SKU: EQFR71

SKU	CAPACITY	HEIGHT x WIDTH x DEPTH	SET POINT RANGE	AMPERAGE	SHIPPING WT.	DEFROST	CERTIFICATIONS
EQFR71	8 cu.ft.	33" x 44" x 24"	36°F to 46°F	1.4	110 lbs	Frost-Free	UL, UL-C, UL-S
EQFR121	13.1 cu.ft.	33.25" x 53.75" x 28.63"	36°F to 46°F	1.49	145 lbs	Frost-Free	UL, UL-C, UL-S
EQFR151	17.1 cu.ft.	33.25" x 65" x 30.25"	36°F to 46°F	1.4	175 lbs	Frost-Free	UL, UL-C, UL-S
EQFR221	23.8 cu.ft.	32" x 83.25" x 31.5"	36°F to 46°F	3.1	210 lbs	Frost-Free	UL, UL-C, UL-S
EQFF72	8 cu.ft.	33.25" x 44" x 24"	-5°F to 12°F	1.4	110 lbs	Frost-Free	UL, UL-C, UL-S
EQFF122	13.1 cu.ft.	33.25" x 53.75" x 28.63"	-5°F to 12°F	1.49	162 lbs	Frost-Free	UL, UL-C, UL-S
EQFF152	17 cu.ft.	33.25" x 65" x 30.25"	-5°F to 12°F	1.4	175 lbs	Frost-Free	UL, UL-C, UL-S
EQFF222	23.8 cu.ft.	32" x 83.25" x 31.5"	-5°F to 12°F	15	210 lbs	Frost-Free	UL, UL-C, UL-S
SCFM62	7 cu.ft.	32" x 38.5" x 22.13"	-22 to 10.4°F	1.4	93 lbs	Manual	UL, UL-C, UL-S
SCFM92	9 cu.ft.	33" x 44" x 24"	-22 to 10.4°F	5	135 lbs	Manual	UL, UL-C, UL-S
SCFM162	14.1 cu.ft.	33.25" x 53.75" x 28.63"	-22 to 10.4°F	1.49	140 lbs	Manual	UL, UL-C, UL-S
SCFM182	18 cu.ft.	33.25" x 65" x 28.63"	-22 to 10.4°F	1.4	197 lbs	Manual	UL, UL-C, UL-S
SCFM232	24.8 cu.ft.	32" x 83.25" x 31.5"	-10 to 10°F	3.1	230 lbs	Manual	UL, UL-C, UL-S
SCFM252W	26.66 cu.ft.	36.75" x 91.25" x 28.25"	-10 to -26°C	3.8	270 lbs	Manual	UL, UL-C, UL-S

Product Options & Accessories

Virtually all Accucold refrigerators and freezers can be configured to meet the unique needs of any office, laboratory, clinic, hospital, pharmacy, or other institution.

NIST Calibrated Temperature Monitoring

Our ISO/IEC 17025:2017 certified laboratory allows us to calibrate temperature monitoring equipment to NIST standards in our own facilities, ensuring accurate temperature readings for your records.



OPTION: DL2BKIT

USB Digital Data Logger

Available installed on your equipment or for aftermarket use, the DL2B DDL with alarm provides virtually unlimited logging memory to record unit temperature to your user-set logging intervals. To help support meeting CDC/VFC guidelines, this logger also keeps the current and high/low temperature in view



OPTION: DLSA11CAL-24P

Portable Loggers

Available in 12 or 24 packs, the DLSA loggers are ideal for cold chain management during the transportation process, with an integrated USB to ensure easy record keeping with downloadable software


New Solutions to Reduce the Spread of Germs



OPTION: STO

Step-To-Open Door Pedal

Factory-installed on select refrigerators and freezers, allowing staff to open a unit's door with a light step



OPTION: TBC

OPTION: CP

Antimicrobial Copper Handles

Towel bar or pipe style handles in pure antimicrobial copper installed on select refrigerators and freezers

Additional Options:



OPTION: DR

Interior Options


- Full-extension ventilated drawers**
- Locking compartments** for individual access
- Access ports** with dry contacts
- 3/8" or 1/2" probe holes**



OPTION: STACKING RACK

Exterior Options

- Temperature alarms** with external readouts
- Self-closing door swings**
- Dollies or NSF-listed 3"/6" locking casters**
- Hospital grade cord** with 'green dot' plugs
- Stainless steel exterior finishing**
- Fan-assist options** for compact units
- Stacking racks** to combine units



OPTION: ALARM

Additional Services

- Temperature displays** in °F or °C
- Door reversals**
- Pyxis-ready doors**



OPTION: PF1KIT

Audible Alarms

- Power-out alarms**
- High/low temperature alarms**
- Open door alarms**



OPTION: KAL

Locking Options

- Keyless Access Lock**
- Battery or AC powered keypad locks**
- Combination lock** with key override



All NIST calibrated equipment

can be re-calibrated after its certification expires by sending the device or probe to our ISO/IEC 17025:2015 certified laboratory

Additional options are available

Contact our product specialists for more information on how we can customize any unit to your facility's precise storage needs

accucold

Learn more online at accucold.com.

For additional assistance, please refer to the directory below.

Sales Assistance:

888.4.MEDLAB ext. 509

info@accucold.com

Service Help:

888.4.MEDLAB ext. 508

service@summitappliance.com

Technical Support:

888.4.MEDLAB ext. 513

techsupport@summitappliance.com

Replacement Parts:

888.4.MEDLAB ext. 512

parts@summitappliance.com

**Medical,
Laboratory,
Pharmacy,
& Life Science
Cold Storage**